

PPI

Quality. Packaging. Solutions.

Gewerbestraße 3 | D-91560 Heilsbronn

Phone +49 (0) 9872 9788-0 | Fax +49 (0) 9872 9788-77 | www.ppi-pharma.de | contact@ppi-pharma.de



Ampoules


Ampullen






Ampoules

Ampollas

Ampoules

Ampullen
Ampoules
Ampollas



 ml	 \pm mm	 h mm Form B	 h mm Form C	 h mm Form D
1	10,75	60	67	70
2	10,75	72	79	83
3	12,75	75	82	89
5	14,75	83	90	95
10	17,75	102	109	112
20	22,5	113	120	126
25	22,5	128	135	141
30	22,5	143	150	156